



Figure 5: CTOR Intrusion-Plates (I-Plates) can be used at different stages of treatment. This CTOR Plate design is used for bilateral intrusion of posterior teeth in the upper arch. The I-Plate (A) can be easily shaped by hand or pliers to follow the patient's palatal contour. This plate can be used in combination with buccal TADs, depending on the patient's problem (B, intrusion forces shown as black arrows). Buttons on the lingual surface of the posterior teeth can be connected to the I-Plate using power thread or power chain (C). The design allows intrusion of several posterior teeth simultaneously. After intrusion this plate can be used for posterior anchorages during retraction of anterior teeth.